

TO Validation

The Tempus Tumor Origin (TO) test uses information from analysis of nucleic acids by next-generation sequencing (NGS) performed as part of a separately-ordered Tempus xR RNA test. The TO test uses tumor mRNA expression results to predict the highest probability primary and other potentially likely cancer types from 68 possible diagnostic cancer types. The Tempus TO test can be reported only for samples with $\geq 20\%$ tumor purity.

Validation of the Tempus TO test was performed in a CLIA-certified, CAP-accredited lab. An analytical validation cohort consisting of 9,210 tumor samples of known origin and 1,708 cancers of unknown primary (CUPs) was created from the Tempus database. The validation cohort consisted of 25% of the labeled samples, which were selected via stratified random sampling, from within the Tempus database. All CUPs were used for the biomarker concordance analysis. All samples with an mRNA expression profile passing the quality control checks were eligible for the study. These samples were assigned one of 68 diagnostic labels or designated as CUPs based on the diagnosis assigned to the sample by Tempus pathologists at the time of sample accessioning and histologic review. Samples included formalin-fixed, paraffin-embedded (FFPE) slides, FFPE tissue blocks, blood, bone marrow aspirates, and fresh frozen tissue.

Labeled samples (n= 9,210) were used to characterize the performance of the TO test. The samples contained a broad range of tumor purities as determined by manual pathologist review, ranging from $< 20\%$ to 100% tumor purity. Performance metrics of the assay were stratified by type and subtype; the overall classification accuracy of the TO model was 91%.

Samples labeled as a CUPs (n= 1,708) were not used in determining performance specifications of the classifier, but were utilized in other studies to characterize reportable range and interassay reproducibility.

The sixty eight (68) possible diagnostic types of the TO test are:

Acute lymphoblastic leukemia	Goblet cell adenocarcinoma	Prostate neuroendocrine carcinoma
Acute myeloid leukemia	Gynecological clear cell carcinoma	Prostatic adenocarcinoma
Adenoid cystic carcinoma	Head and neck squamous cell carcinoma	Renal chromophobe carcinoma
Adrenal cortical carcinoma	Hepatocellular carcinoma	Renal clear cell carcinoma
Anogenital squamous cell carcinoma	High grade glioma	Renal papillary carcinoma
B cell lymphoma	Leiomyosarcoma	Rhabdomyosarcoma
Breast carcinoma	Liposarcoma	Salivary carcinoma
Carcinosarcoma	Low grade glioma	Schwannoma
Cervical carcinoma	Lung adenocarcinoma	Skin neuroendocrine carcinoma
Cholangiocarcinoma	Lung squamous cell carcinoma	Skin squamous and basal cell carcinoma
Chondrosarcoma	Medulloblastoma	Small bowel adenocarcinoma
Chronic lymphocytic leukemia	Melanoma	Small cell lung carcinoma
Chronic myeloid leukemia	Meningioma	Synovial sarcoma
Colorectal adenocarcinoma	Mesothelioma	T cell lymphoma
Endometrial serous carcinoma	Metaplastic breast carcinoma	Thymic squamous cell carcinoma
Endometrial stromal sarcoma	Multiple myeloma	Thyroid cancers
Endometrioid carcinoma	Neuroendocrine lung tumor	Urothelial carcinoma
Ependymoma	Oligodendroglioma	Urothelial neuroendocrine carcinoma
Ewing sarcoma	Osteosarcoma	Vascular sarcoma
Fibrous sarcoma	Ovarian mucinous adenocarcinoma	Well differentiated gastrointestinal neuroendocrine tumor
Gastroesophageal adenocarcinoma	Ovarian serous carcinoma	
Gastroesophageal squamous cell carcinoma	Pancreatic adenocarcinoma	
Gastrointestinal neuroendocrine carcinoma	Pancreatic neuroendocrine tumor	
Gastrointestinal stromal tumor	Peripheral nerve sheath tumor	